



Sterile Fluid Handling Bags 5 Liter, 10 Liter, 20 Liter



TC Tech 10 Liter Standard



TC Tech 20 Liter SCD®

5

TC Tech 5 Liter, 10 Liter, 20 Liter Sterile Fluid Handling Bags accommodate fluid volumes commonly associated with bench-top and pilot plant work. Typical applications include buffer storage and harvest collection. Hanger Ports facilitate complete drainage. May be frozen to -30°C . Consult TC Tech for recommendations on handling frozen bags.

Specifications

TC Tech Standard

- ULDPE Fluid Contact Surface
- Septum Injection Port
- 1 line C-Flex Formulation 082 – 6" Length – 1/4" ID
- 1 line C-Flex Formulation 082 – 6" Length – 3/8" ID
- 2 QDC Coupling Inserts
- 2 Tube Ratchet Clamps
- 2 Port Covers
- Provided Sterile via 25kGy – 40kGy Gamma
- Individually Packaged

TC Tech SCD®

- ULDPE Fluid Contact Surface
- Septum Injection Port
- 1 line C-Flex Formulation 082 – 3' Length 1/8" ID
- 1 line C-Flex Formulation 082 – 6" Length – 3/8" ID
- 1 QDC Coupling Insert
- 1 Male Luer with Female Luer Plug
- 2 Tube Ratchet Clamps
- 1 Port Cover
- Provided Sterile via 25kGy – 40kGy Gamma
- Individually Packaged

Options

Fittings

- QDC Sealing Cap
- QDC Coupling Body / Sealing Plug
- Sanitary Fitting
- Male / Female Luer with Plug
- Hose Barb
- Female Luer for Syringe Connection

Tubing

- Platinum Cured Silicone
- Peroxide Cured Silicone
- Pharmed®
- PVC
- C-Flex® Formulation 001, 050, 072

Other

- Sterilizing Filter
- Manifold System
- Steam-in-Place Bag Connector

**Call for
Custom Design
Consultation**

Ordering Information

Item Number	Item Description	Quantity
111320	Storage Bag, 5 Liter Standard	20 per Case
111420	Storage Bag, 5 Liter SCD	20 per Case
112320	Storage Bag, 10 Liter Standard	20 per Case
112420	Storage Bag, 10 Liter SCD	20 per Case
113315	Storage Bag, 20 Liter Standard	15 per Case
113415	Storage Bag, 20 Liter SCD	15 per Case

